PM Laborhaube AV

Lab Hood AV

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Instruction Manual

1. General Information

The PM Lab Hood AV is designed to protect the breathing organs and the face from gaseous and particulate pollutants. Together with a battery-operated blower motor and filters, the hood can be used as a respiratory protection system in accordance with DIN EN 12941

The head piece can only be operated with an approved air source, as this is the only way to create the required overpressure / minimum nominal volume flow of 120 l/min, which prevents pollutants from entering the PM Lab Hood AV. The air source can be a battery-operated blower filter unit according to EN 12941 or a compressed air control valve according to EN 14594. With an approved blower filter device, the nominal operating time with the PM Lab Hood AV is at least 4 hours.

If you are unsure about the selection and handling of the equipment, please contact your supervisor or retailer. You can also get in contact with the technical customer service of PM Atemschutz, e-breathe Safety or selected service centers.

Areas of use:

The PM Lab Hood AV can be used in combination with a blower filter unit as an alternative to filter protection wherever filtering respiratory protection is recommended.

The PM Lab Hood AV can be used in combination with a compressed air control valve as an isolating device wherever isolating respiratory protection is recommended.

When selecting the appropriate respiratory protection (filtering respiratory protection or insulating respiratory protection), the following factors must be considered:

- Type of Contamination
- Concentration
- Work load
- Need for protection that goes beyond the need for respiratory protection

The risk analysis must be carried out by persons with the relevant training and experience in the area concerned.

1.1 Selection of filter devices with blower

See also DGUV 112-190 (formerly BGR 190 (ZH 1/701)) Rules for the use of breathing apparatus. See also BGI 504-26 (ZH1/600.26) Respiratory protection suitability (medical check-up). See the operating instructions for the individual blower unit and filter.

1.2 Material Characteristics

Version: Limited-Use

The PM AV Lab Hood is made of Polysafe material. For more information on the characteristics of this material, see the Datasheet.

Version: Premium (Reusable):

The PM AV Lab Hood is made of Polysafe material. For more information on the characteristics of this material, see the Datasheet.

1.3 Download / Declaration of Conformity

The currently valid documents (declaration of conformity, certificates, instructions for use, data sheet and brochure) for this product can be found on our homepage under the following link www.e-breathe.de or by sending an e-mail to the following address: info@e-breathe.de

1.4 Information about Certification- & Supervisory Authorities

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2. Restrictions / Terms of Use

GENERAL EN12941 / EN14594

2.1	This product may only be used by qualified persons who are fully aware of the hazards in the workplace and are familiar with the use and handling of the device.
2.2	Respiratory protective equipment must be selected according to the type and concentration of the hazardous substances.
2.3	Insufficient warning properties of hazardous substances require special rules of engagement.
2.4	The tightness requirements cannot be achieved if there is hair or parts of clothing between the hood collar and the neck. The degree of protection is only achieved if the hood is put on correctly!
2.5	When using the respirator in rooms with an explosion hazard, the relevant regulations must be observed.
2.6	The PM Lab Hood AV must be used with the blower unit running. If the blower fails, no or only a low respiratory protection function of the device can be expected. There is a risk of immediate, severe accumulation of carbon dioxide. Blower failure is considered an exceptional condition.
2.7	Ensure that the breathing hose does not get caught on objects in the environment to avoid the risk of the hose or breathing connection becoming detached from the head.
2.8	During heavy breathing under high workloads, negative pressure can develop in the breathing connection during the inhalation phase, which can lead to the penetration of contaminants.
2.9	The protection factor may be reduced if the equipment is used in environments with high wind speeds.
2.10	The PM Lab Hood AV must not be used in case of difficult breathing, if you perceive odor or taste of impurities or you experience any kind of discomfort (dizziness, nausea or similar).
2.11	If you are unsure about the selection and handling of the equipment, please contact your supervisor or the point of sale.

ISOLATING RESPIRATORY PROTECTION EN14594

2.12	The breathing air to be used must comply with EN12021.
2.13	The water content of the breathing air must be kept within the limits of EN 12021 to avoid freezing of the device.
2.14	The product must not be used in atmospheres containing oxygen or oxygen-enriched air.
2.15	The device may not provide adequate protection in highly toxic ambient atmospheres.
2.16	When using hearing protectors or sound-reducing communication equipment, the reduced audibility of the compressed air control valve warning must be taken into account and appropriate instructions must be prepared.
2.17	Minimum volume flow compressed air control valve >160 l/min Maximum volume flow compressed air control valve < 310 l/min

FILTERING RESPIRATORY PROTECTION EN12941

- 2.18 Filter devices should not be used if the environmental conditions are unknown. In cases of doubt, isolating systems should be used which operate independently of the environmental atmosphere.
- 2.19 Filter units are not allowed to be used in confined spaces (e.g. containers, tanks), as a lack of oxygen or heavy oxygen displacing gases (e.g. carbon dioxide) may occur.
- 2.20 When a filter unit is used, the oxygen content of the surrounding air must be at least 17% by volume (this limit may vary according to national regulations). The product cannot be used if the ambient air does not have a normal oxygen concentration or in atmospheres with oxygen or oxygen-enriched air.
- 2.21 A gas filter does not protect against particles and a particle filter does not protect against gases and vapours. If in doubt, use combination filters.
- 2.22 Normal filter devices don't offer protection against certain gases such as CO (carbon monoxide), CO₂ (carbon dioxide) or N₂ (nitrogen) and nitrogen oxides.
- 2.23 Particle filters against radioactive substances and microorganisms (viruses, bacteria, fungi and their spores) should only be used once.
- 2.24 Carefully check the identification of the filters to be used together with the blower and make sure not to confuse the classification according to EN 12941 or EN 12942 with the classification of the blower according to other classification standards.
- 2.25 Filters must not be attached directly to the face shield or breathing hose.
- 2.26 When using the respirator in areas with explosion hazards, the relevant regulations must be observed. The use of the PM Lab Hood AV is only permitted there in conjunction with an explosion-proof (EX) blower unit.
- 2.27 **Only the PM AV Lab Hood with the blower filter unit PM Proflow EX 120 /160 may be used in an EX area!**
- 2.28 **The PM AV Lab Hood with the PM Proflow SC 120/160, e-breathe e-Flow, e-breathe Smartblower blower filter units CANNOT be used in an EX Area!**
- 2.29 Minimum volume flow of blower filter unit >120 l/min
- 2.29 Minimum volume flow PAPR:
- PM Proflow 2 SC 120 > 120 l/min
 - PM Proflow 2 Ex 120 > 120 l/min
 - e-breathe e-Flow: > 160 l/min
 - Smartblower: > 135 l/min
- 2.30 Battery Information PAPR:
- PM Proflow 2 SC 120: 9,6 V / 4.0 AH
 - PM Proflow 2 Ex 120: 9,6 V / 3,8 AH
 - e-breathe e-Flow: 14,4V, 3,4Ah, 49Wh
 - Smartblower: 11,25V, 2,95Ah, Lithium-ion Battery

3. Commissioning

Note: For commissioning, also use the enclosed operating instructions for the compressed air control valve, the blower filter unit and the respiratory protection filters. Only approved components and filters may be used.

3.1 Check before use

Make sure that the equipment is complete and has not been damaged during transport.

- General visual inspection:
 - no cuts, leaks or defects in the material (visor, exhalation valve, hood material) should be visible.
- Check the condition and integrity of the breathing air hose, blower filter device and filters or compressed air control valve and compressed air hose.
- Assemble the breathing air hose to the air source.
- Switch on the air source

Optional Checks:

- **Check the minimum nominal volume flow of the blower** unit according to 2.20. Attach the test tube to the fully charged blower equipped with new filters and switch it on. Read off the air flow, if it is too low, do not use the blower unit.
- **Checking the warning device of the blower:** Close the air outlet of the fully charged blower with the palm of your hand. Shortly thereafter you will notice that the blower is running more intensively or the speed of the blower unit is increasing or readjusting. After some time an alarm should start. If the blower does not give an alarm, the blower should not be used

3.2 Application of the Hood

1. After unpacking, connect the Velcro fastener on the visor with the Velcro fastener on the headband. **Fig. 1**
2. Screw the hose to the air source and apply it according to the instructions for use.
3. Put on the Hood and adjust the headband fixation correctly. The visor should be placed in front of the face. **Fig. 2**
Please note that the outer headband of the hood should be attached under the hose!
4. Attach the air hose to the hood using the click adapter. **Fig.3**
The air hose must run over the back to the blower unit without being warped.
5. Then pull the collar tight with the help of the sewn-in cord. **Fig. 4**
Make sure that there is no hair or clothing between the hood and the neck seal.
6. Tie the outer straps at the front of the hood into a bow. **Fig.5**

3.2 Putting down the Hood

1. Remove the air hose from the hood.
2. Reach into the hood with your hands (do not touch the contaminated material) and open the collar.
3. Pull the hood from the inside off of the head and dispose it properly.

4. Maintenance

Important: Maintenance and cleaning work may only be carried out by appropriately trained personnel who are well trained in this kind of task.

system component	Work to be performed	Before Use	After Use	Semiannually	Two years	Bei Bedarf
Headtop + Hose	View-tight-and functional testing	x		x		
	Control through equipment carrier	x				
	Cleaning and disinfection*		x		x	x
	Exchange of exhalation valve*			x		x
Filter	Expiry date check	x				
	visual inspection	x				
	Check filter capacity	x				x
Blower Unit	See instructions for use of the Blower Unit! (change filter, check/charge battery, cleaning and disinfection and annual service appointments)					
Compressed Air-Control Valve	See operating instructions of the compressed air control valve! (Change seals, check couplings, check warning point, cleaning and disinfection and annual service appointments).					

*) This maintenance applies only to the reusable hoods.

The maintenance operations that can be performed directly by the user are listed below. These maintenance operations are only possible and necessary for the reusable version of the hood.

Disassembly of the exhalation valve

Hold the exhalation valve with one hand outside of the hood. With the other hand, reach inside the hood and turn the exhalation valve cover counterclockwise. Then remove the exhalation valve cover and washer and pull the rest of the exhalation valve out of the hood with the other hand.

Mounting of the exhalation valve

The exhalation valve consists of three parts: the exhalation valve outlet, a washer and an exhalation valve cover. Take the exhalation valve outlet and insert it with the thread from the outside through the hole in the hood. Make sure that the air outlet is directed downwards. Next, turn the canopy so that you can see the thread of the exhalation valve outlet inside the canopy.

Then place the washer on the thread and screw the exhalation valve cover clockwise onto the thread. Hold the exhalation valve from the outside so that it does not slip.

Finally, check that the exhalation valve is correctly seated after installation and that no leakage is visible.

Important: Only use original products from e-breathe. Do not make any changes to the equipment. The use of non-original parts or modifications to the equipment can reduce the protective function, endanger product approvals and permanently damage the equipment. The guarantee expires in case of non-compliance.

5. Cleaning and Disinfection

Note: Cleaning / Disinfection should only be done on the reusable hood.

Important: It is important to be careful during the work steps otherwise the components may be damaged. Therefore, only use the described procedures. Other processes or cleaning agents could damage components.

Important: No solvents (e.g. acetone, turpentine) or bleaching agents (perborate, percarbonate), hot water, compressed air or pressurized water may be used for cleaning.

5.1 Cleaning

Commercial dishwashing detergent can be used for daily care. Dilute the detergent with lukewarm water and wipe off the parts to be cleaned with a cloth. Heavy soiling can be removed carefully with a soft brush. Then wipe dry with a cloth and leave to dry in the air. Protect from direct sunlight.


5.2 Disinfection

The following disinfectants are recommended: PM Desk or Curacid PSA Ultra. Dilute disinfectant with lukewarm water (mixing ratio see imprint bottle) and wipe off the parts to be cleaned with a cloth. Heavy soiling can be removed carefully with a soft brush. Then wipe dry with a cloth and leave to dry in the air. Protect from direct sunlight.

6. Storage

Before storing the reusable hoods, clean the hood in advance. After cleaning, keep the equipment in a dry and clean place, protected from direct heat and sunlight. The blower should be stored in a closed container or cabinet, protected from dust, light, vapours of chemical agents and, if possible, away from heat sources. An appropriately stored, unused Breathing connection is also functional after long storage.

The following storage temperatures should be followed for the hood.

7. Product Label & Packaging

Refer to page 12 for illustrations of the name plate of the hood and packaging.

Nr.	Description		
1	Product Name	10	Symbol Observe Instructions for Use
2	Part. No	11	Symbol CE Marking
3	Approved Standards + Protection Class	12	Do not wash with water
4	Place of Production	13	Do not tumble dry
5	Manufacturer + Address	14	Do not dry clean
6	Version	15	Do not iron
7	Date of manufacture / Serial Number	16	Low flame retardancy
8	Size	17	Number of pieces per packaging unit
9	Symbol Temperature and Humidity		

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Please note: The classification of blower-assisted respiratory protection devices within the scope of the approval may change due to the replacement of system components (blower, head part, hose, etc.). Therefore, please include the latest published technical data in your decision or contact our specialist staff who will be happy to advise you and assist you in putting together the appropriate equipment combination.

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Service Partner
www.pm-atemschutz.de
info@pm-atemschutz.de

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